

CRYSTALHD

Unmatched Brightness.

3x stronger lens. Engineered for lasting clarity.



A brighter reflection with crystal-clear accuracy - now 3x stronger lens for increased scratch resistance.

Less eye strain. Clearer visualization. More precise, confident dentistry.



**Four Lens Sizes.
Standardized for
every procedure.**

Available in #0, #3, #4,
and #5 mirror sizes



**Four Configurations.
Built for your workflow.**

Includes one-piece,
resin, stainless steel,
& double-sided



**16 Colors.
Clear organization.**

Available in multiple
pack sizes.

Trusted by clinicians who prioritize clarity and efficiency.



Contact

Sales: Nichol Schreiber | nichol@zirc.com
INTL Sales: Renata Elez | renata@zirc.com

zirc.com

Crystal HD® Features & Benefits

Unmatched Brightness.

Enhanced Lens Brightness

Delivers exceptional brightness and clarity for improved clinical visualization during exams and procedures.



Improved Scratch Resistance

The advanced lens improves initial surface scratch resistance by 269% compared to previous models, helping preserve the clarity over time.

Standard Rhodium Lens



Crystal HD® Lens



Flush-Mounted Resin Frame

Flush lens integration minimizes buildup and corrosion while eliminating galvanic shock associated with traditional metal mirrors.

”

As a microscope user, I understand the importance of clear vision. Zirc mirrors significantly improved the brightness of my image. Brightness makes all the difference."

Dr. Dimitris Tsanaktsidis

Microscopic & Aesthetic Dentistry

Contact

Sales: Nichol Schreiber | nichol@zirc.com
INTL Sales: Renata Elez | renata@zirc.com

zirc.com



Crystal HD® Features & Benefits cont'd



Precision & Accuracy

Our Crystal HD® Mirror will reduce eye strain with its much brighter and clearer lens reflecting color 99% accurately compared to the 75% accuracy of standard Rhodium mirrors.

- Autoclavable
- No Galvanic Shock
- Flush Lens Design
- Rust-free



Durability & Resistance

Zirc's advanced resin material formula for the Crystal HD® Mirrors provides multiple benefits, including chemical resistance, durability, and patient safety. Resin construction eliminates galvanic shock, a common discomfort associated with traditional metal mirrors, significantly enhancing patient comfort.



Comfort & Ergonomics

The Soft-Grip Handle and the Thin-Grip Handles weighs only **2.3g** and **1.7g** respectively, reducing pinch force and hand fatigue during extensive use.

These handles allow for effortless maneuverability without compromising precise control. Zirc's mirror handles are designed with elastomer grips, offering exceptional comfort and secure handling to minimize strain.

Contact

Sales: Nichol Schreiber | nichol@zirc.com
INTL Sales: Renata Elez | renata@zirc.com

zirc.com