

SDI

HIGHLIGHTS

September - October 2026

50+ YEARS
DENTISTRY
SCIENCE

POLA

ADVANCED TOOTH
WHITENING SYSTEMMADE IN
AUSTRALIA

TAKE-HOME WHITENING

Developed in our state of the art research and development facility and manufactured in Australia, the Pola teeth whitening system continues to set the industry benchmark in speed and performance around the world.

Pola teeth whitening is formulated to safely remove hard, longer term stains from your teeth and return your smile to your whitest natural shade possible. All without damaging the health of your teeth.

Pola Day 6%



POLA 4 X 1.3g MINI KIT	MRP \$41.47
POLA 10 X 1.3g SYRINGE KIT	MRP \$101.57
POLA 50 X 3g DISPENSER PACK	MRP \$595.34
POLA FOR ALIGNERS	MRP \$87.10
POLA LIGHT	MRP \$114.61
POLA LUMINATE	MRP \$28.80

BUY 4
GET 1
FREE^x



IN-OFFICE WHITENING



SUPER. FAST.
WHITENING.

POLA RAPID

Pola Rapid is an in office whitening treatment that safely whitens teeth in a 24 minute application time.

SUPERIOR HANDLING AND FAST REMOVAL

Pola Rapid's higher viscosity gel enables faster application and removal due to its specially formulated non-stick attributes. Pola Rapid can be extruded quickly and easily and remains on teeth when applied.

ENHANCED BLUE GEL FOR FASTER APPLICATION

Pola Rapid uses an advanced blue gel for ease of visibility and application. The easy to apply non-stick gel enables faster application and removal.

MINIMAL CHAIR TIME WITH IMMEDIATE RESULTS

A unique **38% hydrogen peroxide** formulation that releases active ingredients immediately. In only one visit, Pola Rapid initiates a powerful whitening treatment designed to whiten teeth in 24 minutes. No pre-desensitiser step is required.

BUY A POLA RAPID KIT
(1 PATIENT OR 3 PATIENT KIT)

**GET A BONUS
POLA LIGHT**
FOR EACH PATIENT



POLA RAPID 3 PATIENT KIT	MRP \$473.78
POLA RAPID 3 PATIENT KIT WITH OPTRAGATE	MRP \$504.63
POLA RAPID 1 PATIENT KIT	MRP \$167.93
POLA RAPID 1 PATIENT KIT WITH OPTRAGATE	MRP \$179.67



STELA

THE FUTURE OF COMPOSITES

50 years of dental material experience has engineered a two step restorative that is stronger, gap-free, easier to place and perfectly cures every time.

- ✓ Unlimited depth of cure (self cure)
- ✓ Gap-free interface
- ✓ Only 2 steps: 15 sec prep
- ✓ High strength
- ✓ Choose: capsule or automix

STELA CAPSULE REFILL 50 CAPS	MRP \$291.98
STELA CAPSULES INTRO KIT	MRP \$145.54
STELA AUTOMIX REFILL	MRP \$200.82
STELA AUTOMIX INTRO KIT	MRP \$243.24

**BPA &
HEMA
FREE**



LEARN MORE

25% OFF

STELA AUTOMIX

ADHESIVES

ZIPBOND

THE ALL-IN-ONE ADHESIVE SOLUTION FOR EVERY SITUATION

UNIVERSAL ADHESIVE

- Consistent bond strengths to dentine and enamel
- Wet and dry dentine bonding
- Proprietary combination of adhesive monomers including MDP
- Predictable bond strengths to indirect substrates
- Effective hybridisation
- Exceptional control
- Highly aesthetic
- Pleasant aroma

ZIPBOND REFILL 5mL BOTTLE MRP \$162.95

25% OFF*



29 CLINICAL EVALUATORS
862 TOTAL USES
93% CLINICAL RATING

MDP
MONOMER

BPA &
HEMA
FREE

ACCESSORIES

POINTS


DISPOSABLE BRUSH APPLICATORS

- Points disposable brush applicators are used to quickly transfer small amounts of material into areas that are difficult to reach. Points are resistant to common solvents in dental materials.
- Solvent resistant, non-absorbent, non-linting fibers
- Boxes of 400 in medium, fine and superfine
- The three sizes of Points, superfine, fine and medium, enable one to pick up different drop sizes.
- The various colors are ideal for color coding multi-step procedures.

BOX OF 400 POINTS MRP \$51.91
POINTS DISPENSER MRP \$43.07



Points can be used either
straight or bent



BUY 2, GET 1 FREE*

THE NEXT GENERATION SDF SYSTEM

SDF INNOVATION: RIVA STAR

RIVA STAR

Drawing on 50 years of dental research experience and proudly developed in SDI's state of the art manufacturing facility in Australia, Riva Star is a silver diamine fluoride solution system used to desensitise tooth pain immediately.

This next generation SDF treatment delivers Minimal Intervention Dentistry that adopts the philosophy of integrating prevention, remineralisation and minimal intervention.

RIVA STAR BOTTLE KIT	MRP \$144.36
RIVA STAR CAPSULE KIT	MRP \$131.23



LEARN MORE

NEXT GENERATION SDF: RIVA STAR AQUA

RIVA STAR AQUA

SDI continues to innovate SDF solutions by releasing Riva Star Aqua, a next generation water-based silver fluoride solution that is also effective as an SDF solution, with additional patient benefits.

Similar to Riva Star (SDF), Riva Star Aqua [AgF] is also a non-invasive 2 step patented system used to desensitise tooth pain in a few drops. Riva Star Aqua is an aqueous silver fluoride solution without the Ammonia base.

RIVA STAR AQUA BOTTLE KIT	MRP \$158.79
RIVA STAR AQUA CAPSULE KIT	MRP \$144.36



LEARN MORE



\$110^{*}

BOTTLE KITS

\$125^{*}

CAPSULE KITS



ADVANTAGES	RIVA STAR	RIVA STAR AQUA
Two years and beyond desensitising effect	✓	✓
Staining reduction / elimination	✓	✓
Deep penetration and prevention of biofilm	✓	✓
Antimicrobial activity	✓	✓
Increasing bond strengths of GICs	✓	✓
Storage requirements	Bottles: Room temperature Capsules: Refrigerated	✓ Room temperature
Soft tissue irritation	Mild discomfort	✓ No discomfort
Odour / taste	Slight ammonia taste and odour	✓ No ammonia taste or odour
Non-invasive / MID	✓	✓

¹ Turton B, Horn R, Durward C. Caries arrest and lesion appearance using two different silver fluoride therapies on primary teeth with and without potassium iodide : 12 month results. Clin Exp Dent res 2020;1-11

EQUIPMENT

ULTRAMAT 2 \$750*

The Ultramat 2 offers extremely precise, consistent and accurate trituration. Furthermore, it can be used with most encapsulated dental materials. With a high energy output, it has a mixing frequency of 4600 oscillations per minute but uses remarkably little power. Best of all due to its robust and reliable construction, the 5 year warranty assures your purchase is a long lasting investment.

The multi-use Ultramat 2 is suitable for all types and brands of glass ionomer, amalgam and cement capsules.

The micro-switch automatically halts trituration when the lid of the Ultramat 2 is lifted.



SAVE
\$862.70*

ULTRAMAT 2 MRP \$1612.70

5 YEAR
WARRANTY

1-16
SECOND
MICROPROCESSOR

4600
OSCILLATIONS

**120 &
240V~**

PLASTVAC \$950*

TRAYFORM VACUUM MACHINE

- Powerful motor
- Quick and practical
- Suitable for round and square thermoplastic materials of different dimensions

PLASTVAC MRP \$1614.60



SAVE
\$664.60*

HENRY SCHEIN®

Rely on Us™

orders phone 1300 65 88 22 | orders fax 1300 65 88 10

SDI

YOUR
SMILE.
OUR
VISION.

Call 1800 337 003
info@sdi.com.au

www.sdi.com.au
www.polawhite.com.au

New South Wales
Elli Mouhtouris 0400 849 017
Aurelie Lequimener 0418 210 341
Susanna Kitchen 0400 305 512

Western Australia
Toni Doherty 0448 086 925
South Australia/Northern Territory
Jade Knowles 0417 080 771

Victoria
Jo Coyle 0419 313 765
Elly Cheraghvandi 0418 567 758
Holly Strampel 0408 164 626

Queensland
Natalie Johnson 0409 807 290
Dora Zanin 0428 910 805
Tasmania
Elly Cheraghvandi 0418 567 758

*Offers are valid 1st September - 31st October 2026. Offers cannot be combined. All prices include GST. Free goods will be of equal or lesser value. Purchases must be on one order to qualify for free goods. The above-mentioned offers are not applicable on corporate, government accounts or those customers who have received discounted pricing. Offers are valid in Australia only.